Xpert Bladder Cancer Panel C104

INTENDED USE:

Xpert Bladder Cancer Panel C104 is intended for use as a quality control to monitor the performance of *in vitro* laboratory nucleic acid procedures for the qualitative detection of five RNA transcripts, CRH, IGF2, UPK1B, ANXA10, ABL1, performed on Cepheid GeneXpert® Systems with assays Xpert® Bladder Cancer Monitor and Xpert® Bladder Cancer Detection.

The American Cancer Society estimates over 70,000 new cases of bladder cancer for 2016 and it is currently the fourth most common cancer among men in the United States¹. Bladder cancer also has the highest reoccurrence rate of any malignancy; therefore, diligent and accurate detection and monitoring is crucial for early diagnosis and treatment.

PRODUCT SUMMARY and PRINCIPLE:

The Xpert Bladder Cancer Panel C104 kit contains Xpert Bladder Cancer Low Positive Control and Xpert Bladder Cancer Negative Control. Xpert Bladder Cancer Low Positive Control contains synthetic RNA transcripts carrying target sequences of CRH, IGF2, UPK1B, ANXA10 and ABL1 designed to be detected by Xpert Bladder Cancer assays. Xpert Bladder Cancer Negative Control contains synthetic RNA transcripts carrying target sequence of the endogenous control ABL1 only.

Routine use of quality controls that are consistent lot to lot assists the laboratory in identifying shifts, trends, and increased frequency of random errors caused by variations in the test system, such as failing reagents. Early investigation can prevent failed assay runs.

COMPOSITION:

Xpert Bladder Cancer Panel C104 contains 6 bottles of Xpert Bladder Cancer Low Positive Control and 6 bottles of Xpert Bladder Cancer Negative Control, 4mL each, composed of synthetic RNA transcripts suspended in a solution of synthetic urine, buffers, preservatives, and stabilizers. Table 1 lists the expected results for the Xpert Bladder Cancer Panel C104 when tested by Xpert Bladder Cancer assays on the GeneXpert Systems.

STORAGE and STABILITY:

Xpert Bladder Cancer Panel C104 should be stored at -20°C. Unopened Xpert Bladder Cancer Panel C104 material is stable through the expiration date printed on the kit label when continuously stored frozen. Xpert Bladder Cancer Low Positive Control and Xpert Bladder Cancer Negative Control are for single use. Discard after use according to your local and federal regulations.

PRECAUTIONS and WARNINGS:

Note: This product does NOT contain any biological material of human or animal origin.

- Use control as provided. Do not dilute.
- For Research Use Only. The performance characteristics of this product have not been established.
- Xpert Bladder Cancer Panel C104 contains guanidinium chloride which is toxic if swallowed (H302), irritating to eyes and skin (H315, H319), and is harmful to aquatic life (H402). See Regulation EC 1272/2008 (EU CLP).
- Xpert Bladder Cancer Panel C104 cannot be cloned, sold, or transferred without the explicit written consent of MMQCI.

¹ American Cancer Society. Cancer Facts & Figures 2016. Atlanta, Ga: American Cancer Society; 2016.

INSTRUCTIONS FOR USE:

1. Allow the control to be tested to come to room temperature (18° – 25°C), approximately 30 minutes. Control should feel warm, not cool, to the touch prior to use.
2. Use the control as provided. **DO NOT DILUTE.**
3. Immediately before use, mix the control by inverting at least 10 times. Tap the bottle 3 times on the bench to remove any control caught in the cap before opening the bottle.
4. Follow the assay manufacturer's Instructions for Use to process the control as a clinical sample.
5. Discard control bottle after use according to your local and federal regulations.

EXPECTED VALUES:

Xpert Bladder Cancer Panel C104 does not have assigned values. It is recommended that each laboratory establish its own performance guidelines for this product.

Table 1: Xpert Bladder Cancer Panel C104 Test Result

Kit Component:	Test Result:
Xpert Bladder Cancer Low Positive Control	Positive
Xpert Bladder Cancer Negative Control	Negative

ORDERING INFORMATION

Xpert Bladder Cancer Panel C104

Part Number: C104

Kit Contains: 12 bottles x 4mL
6 each of Low Positive & Negative Controls

(This page intentionally left blank)

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