



INTROL™ PGx 1 Control Warfarin Genotype Table

Polymorphism (Allele)	Bottle b	Bottle c	Bottle d
CYP2C9			
430C>T (*2)	C/T	C/C	T/T
1075A>C (*3)	A/C	C/C	A/A
1080C>G (*5)	C/G	C/C	G/G
818delA (*6)	A/delA	A/A	delA/delA
1003C>T (*11)	C/T	C/C	T/T
374G>A (*14)	G/A	A/A	G/G
485C>A (*15)	C/A	C/C	A/A
895A>G (*16)	A/G	A/A	G/G
VKORC1 Promoter			
-1639G>A	G/A	G/G	A/A
CYP4F2			
1347G>A	G/A	G/G	A/A