



**INTROL™ PGx 1 Control Warfarin Genotype Table**

<b>Polymorphism (Allele)</b>	<b>Bottle b</b>	<b>Bottle c</b>	<b>Bottle d</b>
<b>CYP2C9</b>			
<b>430C&gt;T (*2)</b>	<b>C/T</b>	<b>C/C</b>	<b>T/T</b>
<b>1075A&gt;C (*3)</b>	<b>A/C</b>	<b>C/C</b>	<b>A/A</b>
<b>1080C&gt;G (*5)</b>	<b>C/G</b>	<b>C/C</b>	<b>G/G</b>
<b>818delA (*6)</b>	<b>A/delA</b>	<b>A/A</b>	<b>delA/delA</b>
<b>1003C&gt;T (*11)</b>	<b>C/T</b>	<b>C/C</b>	<b>T/T</b>
<b>374G&gt;A (*14)</b>	<b>G/A</b>	<b>A/A</b>	<b>G/G</b>
<b>485C&gt;A (*15)</b>	<b>C/A</b>	<b>C/C</b>	<b>A/A</b>
<b>895A&gt;G (*16)</b>	<b>A/G</b>	<b>A/A</b>	<b>G/G</b>
<b>VKORC1 Promoter</b>			
<b>-1639G&gt;A</b>	<b>G/A</b>	<b>G/G</b>	<b>A/A</b>
<b>CYP4F2</b>			
<b>1347G&gt;A</b>	<b>G/A</b>	<b>G/G</b>	<b>A/A</b>