

FILMARRAY® TF Control Panel

Ensure Accurate Diagnostics



Vector-borne tropical infections, previously unique to tropical regions, are now quickly spreading to other parts of the world. Presenting with similar non-specific symptoms, tropical fevers pose diagnostic challenges and require different targeted treatments. Highly accurate molecular testing results are crucial for timely therapy decisions.

MONITOR THE DETECTION OF ALL 6 TARGETS WITH ONE SINGLE TUBE

- ✓ Chikungunya virus
- ✓ Dengue virus
- ✓ *Plasmodium* spp.
- ✓ *Plasmodium falciparum*
- ✓ *Plasmodium vivax/ovale*
- ✓ *Leptospira* spp.

Kit contains:

- 6 Positive Controls
- 6 Negative Controls
- 200 µL each



Reliable:
Lot to lot consistency tested over multiple sites, instruments, operators, and pouch lots



Multiplex:
Monitor the detection of all pathogens included in the assay with a single tube



Ready to use:
Use as provided. No reconstitution, pooling, or dilution needed



Safe to handle:
Does not contain any biological material from human or animal origin



Versatile:
Can be used as a reference material for validation, verification, and routine monitoring