

# Unyvero LRT/Pneumonia Control Panel M404



## RELIABLE

Tested over multiple product lots and pouch lots, at multiple clinical sites with multiple operators and instruments



## SAFE

Composed of synthetic DNA suspended in a non-infectious solution of buffers, preservatives and stabilizers



## MULTIPLEX

Only 2 positive control tubes (+ 1 negative) to monitor the detection of all pathogens and antibiotic resistance markers included in the assay

### Kit contains:

**4 tubes x 180µL of each**

- **Positive A**
- **Positive B**
- **Negative**

## For use with:

 **Unyvero LRT**

 **Unyvero LRT BAL**

 **Unyvero HPN**

Maine Molecular Quality Controls, Inc. 23 Mill Brook Rd Saco, ME 04072